

Question				Marking details	Marks Available					
					AO1	AO2	AO3	Total	Maths	Prac
5				Indicative content: 1. Use the small funnel to fill the burette with acid. Run a little acid out into a waste beaker to fill the part of the burette that is below the tap. Record the starting volume of acid in the burette. 2. Accurately measure 25cm³ of sodium hydroxide solution into a conical flask. 3. Add 2 drops of indicator. 4. Add 0.1cm³ of acid at a time, swirl the flask after each acid addition. Keep adding acid until the indicator changes colour. Record the final volume of acid in the burette. 5. Repeat steps 1-4 twice more. Hydrochloric acid and sodium hydroxide are irritants/corrosive. Wear gloves and wash hands immediately after contact with solutions. Wear goggles to prevent acid splashing into eyes when filling burette. 5 – 6 marks Detailed description of method, hazards and precautions. <i>There is a sustained line of reasoning which is coherent, relevant, substantiated and logically structured. The candidate uses appropriate scientific terminology and accurate spelling, punctuation and grammar.</i>	6			6		6