

Question				Marking details (QER)	Marks available					
					AO1	AO2	AO3	Total	Maths	Prac
6	(a)			Indicative content: - Working close to Bunsen flame/disinfect bench - Label the base of the petri dish(es) to indicate antibiotics - Flame forceps (then cool) - Pick up each antibiotic disc in turn and place on agar surface - Minimum lifting of lid /Seal dish with tape - Incubate {for 2-3 days/ at 20 – 25 °C} (Observe results and) measure diameter of clear area around each disc - compare the results for the antibiotic 5-6 marks At least seven correct points from indicative content <i>There is a sustained line of reasoning which is coherent, relevant, supported by evidence and logically structured. The candidate uses appropriate scientific terminology and accurate spelling, punctuation and grammar.</i> 3-4 marks At least four correct points from indicative content <i>There is a line of reasoning which is partially coherent, largely relevant, supported by some evidence and with some structure. The candidate uses mainly appropriate scientific terminology and some accurate spelling, punctuation and grammar.</i>	6			6		6